

# LIST OF COMPLIANT DISPENSARY REFRIGERATORS



Quality Care  
Pharmacy Program

Supporting Excellence in Pharmacy

## Purpose

QCPP recognises purpose-built dispensary refrigerators on this list as meeting the requirement of 'Compliant Dispensary Refrigerator' in Element 5 Action 4.

Where a pharmacy receives a non-conformance at QCPP Assessment because the dispensary refrigerator(s) is not on this list, we encourage submission of a Remedial Review Request at [www.qcpp.com/resources/tools/qcpp-forms](http://www.qcpp.com/resources/tools/qcpp-forms).

Company	Current compliant models	Discontinued models which are still compliant
<b>All-Lab Scientific</b>		ALV-358
<b>Arcus Australia Pty Ltd</b> 08 9381 3888   <a href="http://www.arcus-au.com.au">www.arcus-au.com.au</a>	SV and PV	
<b>Arrowsmith and Grant Refrigeration</b> 03 9706 6655 <a href="http://www.arrowsmithgrant.com.au">www.arrowsmithgrant.com.au</a>	AG 220 V	
<b>Bromic MediFridge®</b> 1300 276 642 <a href="http://www.bromicrefrigeration.com.au">www.bromicrefrigeration.com.au</a>	MED0140GD, MED0374GD	VSS 374 T
<b>Crown Scientific</b>		<b>Dometic:</b> MP 285, MP 155 & MP 405
<b>CSK Sales (Pinnacle)</b> 07 3376 5277   <a href="mailto:scientific@cskgroup.com.au">scientific@cskgroup.com.au</a> <a href="http://www.cskgroup.com.au">www.cskgroup.com.au</a>	CSK PPR 430, CSK PPR 145, PPR-S-300-1-GD	
<b>Elk Scientific</b> 0423 520 302   <a href="mailto:ask@elkscientific.com.au">ask@elkscientific.com.au</a> <a href="http://www.elkscientific.com.au">www.elkscientific.com.au</a>	<b>Logicold Pharmacy:</b> ELKLCP30, ELKLCP60 ELKLCP125, ELKLCP300 ELKLCP400	
<b>Euro Chill</b> 1300 459 140   <a href="http://www.eurochill.com.au">www.eurochill.com.au</a>	VS80 EC-80 Ltr VS311 EC-310 Ltr, VS 370 EC-370 Ltr, Vacc Safe 80 PLUS, Vacc Safe 250 PREMIUM, Vacc Safe 311 PLUS, Vacc Safe 400 PREMIUM, Vacc Safe 371 PLUS, Vacc Safe 120 PREMIUM	VS400 EC, VS250 EC, VS Select 370 EC, Vacc S Select 311 EC, Vacc Safe Premium 120 EC

Company	Current compliant models	Discontinued models which are still compliant
<b>ICS Pacific</b> 03 9706 5115   <a href="mailto:info@icspacific.com.au">info@icspacific.com.au</a>   <a href="http://www.icspacific.com.au">www.icspacific.com.au</a>	Pharma 1000GD, Pharma 3000GD, Pharma 3000SD	G105L Medical, G125L Medical, G130L Medical, G135L Medical, G315L Medical, G375L Medical, G400L Medical, G410L Medical
<b>LabGear Australia</b> 07 3368 1112 <a href="http://www.labgearaustralia.com.au">www.labgearaustralia.com.au</a>	<b>Pharm Chill:</b> BC-29, BC-140, BC-300, BC-400	
<b>Laboratory Equipment Pty Ltd</b> 02 9560 2811   <a href="http://www.labec.com.au">www.labec.com.au</a>	LABEC PVR-130GD, PVR-374, LABEC PVR - 440GD	
<b>Liebherr Mediline (Andi-Co Australia Pty Ltd)</b> 1 Stamford Road, Oakleigh VIC 3166 <a href="mailto:sales@andico.com.au">sales@andico.com.au</a>   1800 685 899 <a href="http://www.liebherrprofessional.com.au">www.liebherrprofessional.com.au</a>	LKUv1610, LKUv1613, LKv3910, LKv3913	LKUv1612, LKv3912
<b>MediLine Equipment</b> 16 Forum Place, Tullamarine VIC 3043   0400 185 682	125ML, 300ML, 400ML; Labcold: EC110, EC300, EC400	
<b>Nuline Refrigeration</b> 7 Longford Court Springvale VIC 3171 03 9562 3355 <a href="mailto:info@nulinerefrigeration.com.au">info@nulinerefrigeration.com.au</a>	HR200, HR400, HR400P, BEV60P & HR600	TTGPR14L, BEV32P, BEV40P
<b>Orford Refrigeration Pty Ltd</b> PO Box 6085, Toowoomba QLD 4350 07 4639 9000   <a href="http://www.orford.com.au">www.orford.com.au</a>	V V D P07A (small), V V D P09A (med), V V D P14A (large)	
<b>Quantum Scientific</b>		<b>Sanyo Medi-Cool:</b> MPR-161D(H), MPR-311D(H)
<b>Avem Quirks Medisafe Refrigerator</b> 1800 658 120   <a href="http://www.avemquirks.com.au">www.avemquirks.com.au</a>	<b>Medisafe:</b> HB80, 311, 371 (glass or solid door)	<b>Medisafe:</b> R88 (small), S5, 125, 410, 471
<b>Rapini Pty Ltd</b> 03 9555 0848   <a href="http://www.rapini.com.au">www.rapini.com.au</a>	R155-G, R300-G	

Company	Current compliant models	Discontinued models which are still compliant
<b>Rollex Vaccine Refrigerators</b> Rollex Group (Australia) Pty Ltd 1300 880 441 <a href="http://www.rollexmedical.com.au">www.rollexmedical.com.au</a>	<p><b>LEC:</b>            PE109; PE102; PG102;            PE207; PG207; PSR151;            PGR151; PSR273; PGR273;            PSR353; PGR353; PE1607;            PG1607</p> <p><b>Pol Eko - MATOS<sup>1*</sup>:</b> MATOS-eco 68 R/G; MATOS-eco 150 R; MATOS-eco 150 R/G; MATOS-eco 200 R; MATOS-eco 300 R; MATOS-eco 300 R/G; MATOS-eco 493 R; MATOS-eco 493 R/G; MATOS-eco 625 R; MATOS-eco 625 R/G</p> <p><b>Aria – MATOS</b>            MATOS Aria-Eco 220 glass            MATOS Aria-Eco 300 Glass            MATOS Aria-Eco374 Glass</p> <p><b>Angelantoni:</b>            FRL 360; FRL 360 Glass; FRL 500; FRL 500 Glass; 700 TN; 700 TN-GL; 1500 TN; 1500 TN-GL</p> <p><b>Twinbird Portable:</b>            SC-BV25</p>	<p><b>LEC</b>            PE202; PG202; PE 501; PG 501; PE502; PG502; PE 901; PG901; PE902; PG902; PE1202; PG1202; PE1602; PG1602; PE507; PG507; PE907; PG907; PE357; PG357; PSR334; PGR334</p> <p><b>Pol Eko - MATOS:</b>            CHL 1+ A; CHL 2+; CHL 2+ A; CHL 3+ A; CHL, 5+; CHL 5+ A; CHL 500; CHL 500 A; CHL 700; CHL 700 A</p> <p><b>RK - Koldpia:</b>            RK600s; RK600g; RK1100s; RK1100g; SGRR-23SD; SGRR-23GD; SGRR-49SD; SGRR-49GD</p>
<b>Skope Industries</b>		B130-PH, B170E-PH, B310E-PH
<b>Thermoline Scientific Refrigerators</b> 02 9604 3911   <a href="http://www.thermoline.com.au">www.thermoline.com.au</a> <i>Note:</i> 1-GD = 1 glass door 2-GD = 2 glass doors 3-GD = 3 glass doors	TEPR-80-1-GD, TEPR-145-1-GD, TEPR-200-1-GD, TEPR-230-1-GD, TEPR-440-1-GD, TEPR-500-1-GD, TEPR-740-2-GD, TLMPR 220-1-GD	TEPR-140-1-GD
<b>Vaccine Fridges Australia (Frigidum Pty Ltd)</b> 1300 272 464 <a href="http://www.vaccinefridges.com.au">www.vaccinefridges.com.au</a>	Huxford VXH20	
<b>VWR International Pty Ltd</b> 1300 727 696   <a href="mailto:csaus@au.vwr.com">csaus@au.vwr.com</a>	<p><b>Panasonic:</b>            MPR-S163-PE,</p> <p><b>Dometic:</b> MP 285,            MP 155, MP 405</p> <p><b>Sanyo Medi-Cool:</b> MPR-161D (H)            MPR-311D (H)</p>	

<sup>1</sup> \*MATOS-cloud models are a variation of MATOS-eco, which include online monitoring capability and are also compliant.

Company	Current compliant models	Discontinued models which are still compliant
<b>Yili Ramli &amp; Co. Pty Ltd</b> 02 9679 9488 <a href="http://www.vaccinerefrigeration.com.au">www.vaccinerefrigeration.com.au</a>	VaccineCab YC-280	



**Suppliers**

Manufacturers wishing to have a dispensary refrigerator added to this list should contact the Cold Chain Testing Centre on 02 9467 7110. Inclusion on this list does not constitute endorsement or recommendation of a company or their products by QCPP.