

# QCPP Compliant Dispensary Refrigerators



Quality Care  
Pharmacy Program

QCPP recognises purpose-built dispensary refrigerators on this list as meeting the requirement of 'Compliant Dispensary Refrigerator'.

QCPP acknowledges that refrigerators that do not appear on this list may still comply with the QCPP requirements. If a pharmacy receives a Corrective Action at a QCPP Assessment because the dispensary refrigerator(s) is not on this list, we encourage submission of a Corrective Action Review via the QCPP Customer Portal

Company	Compliant Models		
<b>All-Lab Scientific</b>	ALV-358 (discontinued)		
<b>Arcus Australia Pty Ltd</b> 08 9381 3888 <a href="http://www.arcus-au.com.au">www.arcus-au.com.au</a>	SV PV AUD1-SV AUD1-V CMU		
<b>Arrowsmith and Grant Refrigeration</b> 03 9706 6655 <a href="http://www.arrowsmithgrant.com.au">www.arrowsmithgrant.com.au</a>	AG 220 V		
<b>Avem Quirks Medisafe Refrigerator</b> 1800 658 120 <a href="http://www.avemquirks.com.au">www.avemquirks.com.au</a>	<b>Medisafe:</b> Medisafe HB80  <b>Medisafe Plus:</b> Medisafe Plus 80 Medisafe Plus 311 Medisafe Plus 410 Medisafe Plus 371 Medisafe Plus 650 Medisafe Plus 1000  <b>Medisafe Discontinued</b> Medisafe 80 Medisafe 125 Medisafe 311 Medisafe 371 Medisafe 410 Medisafe 471 Medisafe R88 Medisafe S5	<b>AQ Medical</b> AQM 110 AQM 280 AQM 350	<b>Vestfrost</b> BC85 SLC125 AKG1365 SC60
<b>Aurora BioScience Pty Ltd</b> 1300 309 994 <a href="http://www.aurorabioscience.com.au">www.aurorabioscience.com.au</a>	<b>Helmer:</b> iBR105-GX      iLR105-GX      iPR105-GX iBR113-GX      iLR113-GX      iPR113-GX iBR120-GX      iLR120-GX      iPR120-GX iBR125-GX      iLR125-GX      iPR125-GX iBR245-GX      iLR245-GX      iPR245-GX iBR256-GX      iLR256-GX      iPR256-GX HBR105-GX      HLR105-GX      HPR105-GX HBR113-GX      HLR113-GX      HPR113-GX HBR120-GX      HLR120-GX      HPR120-GX HBR125-GX      HLR125-GX      HPR125-GX HBR245-GX      HLR245-GX      HPR245-GX		

	HBR256-GX	HLR256-GX	HPR256-GX
<b>BD Rowa</b> 1800 833 518 <a href="http://www.carefusion.com.au">www.carefusion.com.au</a>	Rowa Vmax system		
<b>Bio-Strategy Pty Limited*</b> 1800 008 453 <a href="http://www.bio-strategy.com">www.bio-strategy.com</a>  <i>*Previously VWR International Pty Ltd</i>	<b>PHCbi*:</b> MPR-161D(H) (discontinued) MPR-311D(H)** MPR-S163-PE*** MPR-514(H)^ MPR-514D(H) MPR-514R^^ MPR-S300H-PE MPR-S150H-PE MPR-S500H-PE MPR-S500RH-PE  <i>*Previously branded Panasonic and Sanyo prior to that. All three brands with above model numbers are approved.</i> <i>**replaced by MPR-S300H-PE</i> <i>***replaced by MPR-S150H-PE</i> <i>^replaced by MPR-S500H-PE</i> <i>^^replaced by MPR-S500RH-PE</i>	<b>Domestic:</b> MP 155 (discontinued) MP 285 (discontinued) MP 405 (discontinued)	
<b>Bromic MediFridge®</b> 1300 276 642 <a href="http://www.bromicrefrigeration.com.au">www.bromicrefrigeration.com.au</a>	MED0140GD MED0374GD MED0220GD		
<b>B&amp;R Refrigeration Solutions</b> 0430 280 865 <a href="http://www.brrefrigerationsolutions.com">www.brrefrigerationsolutions.com</a>	BRSC 125 BRSC 300 BRSC 400 HYC-L130		
<b>Coolsafe Medical</b> <b>1800 101 033</b> <a href="http://www.yesmed.com.au">www.yesmed.com.au</a> <a href="mailto:yes@yesmed.com.au">yes@yesmed.com.au</a>	CS145 CS280 CS350 CS381		
<b>Crown Scientific</b>	<b>Domestic:</b> MP 150 MP 155 (discontinued) MP 285 (discontinued) MP 405 (discontinued)		
<b>CSK Sales (Pinnacle)</b> 07 3376 5277 <a href="http://www.cskgroup.com.au">www.cskgroup.com.au</a>	PR-S-36-1-GD PR-S-66-1-GD PR-S-145-1-GD PR-S-300-1-GD PR-S-430-1-GD PR-M-1081-2-GD  <i>Please note codes updated 9/1/2023, previous codes starting with PPR still acceptable</i>	CSK PPR 145 CSK PPR 430	
<b>Elk Scientific</b> 0423 520 302 <a href="http://www.elkscientific.com.au">www.elkscientific.com.au</a>	<b>Logicold Pharmacy:</b> ELKLCP30 ELKLCP60 ELKLCP125 ELKLCP300 ELKLCP400		
<b>Eco Jemss</b> <a href="https://ecojemss.com.au">https://ecojemss.com.au</a> (domestic site)	EKT80 EKT150 EKT175		

(08) 9304 4361	EKT1450 EKT250 EKT425 EKT725		
<b>Euro Chill</b> 1300 459 140 <a href="http://www.eurochill.com.au">www.eurochill.com.au</a>	Vacc Safe 60/Botox BT60 Vacc Safe 135/Botox BT120 Vacc Safe 80 PLUS Vacc Safe 311 PLUS Vacc Safe 371 PLUS Vacc Safe 120 PREMIUM Vacc Safe 250 PREMIUM Vacc Safe 400 PREMIUM Vacc Safe 600 PREMIUM Vacc Safe 700 PREMIUM Vacc Safe 1300 PREMIUM Vacc Safe 1000 PREMIUM Vacc Safe Premium 120 EC Vacc Safe VS80 EC Vacc Safe VS311 EC Vacc Safe VS371 EC Vacc Safe Select 311 EC Vacc Safe Select 370 EC	VS80 EC-80 Ltr VS311 EC-310 Ltr VS370 EC-370 Ltr VS250 EC VS350 VS350P VS400 EC VS420P	<b>Mediguard:</b> Mediguard 60 Mediguard 110 Mediguard 130 Mediguard 135 Mediguard 315 Mediguard 401 Mediguard 601 Mediguard 1001 Mediguard 300 Premium Mediguard 400 Premium Mediguard 401 Plus Mediguard 601 Plus Mediguard 1001 Plus
<b>Exquisite Marketing Australia</b> 03 9561 1595	MV30 MV60 MV150 MV300 MV400		
<b>ICS Pacific</b> 03 9706 5115 <a href="http://www.icspacific.com.au">www.icspacific.com.au</a>	<b>Discontinued models:</b> G105L Medical G125L Medical G130L Medical G135L Medical G315L Medical G375L Medical G400L Medical G410L Medical	Pharma 1000GD Pharma 3000GD Pharma 3000SD Pharma 3	F135 F375P
<b>LabGear Australia</b> 07 3368 1112 <a href="http://www.labgearaustralia.com.au">www.labgearaustralia.com.au</a>	<b>Pharm Chill:</b> BC-29 BC-140 BC-300 BC-400		
<b>Laboratory Equipment Pty Ltd</b> 02 9560 2811 <a href="http://www.labec.com.au">www.labec.com.au</a>	LABEC PVR-130GD LABEC PVR-440GD LABEC PVR-374 LABEC GMO374 LABEC HR200GP		
<b>Liebherr Mediline (Andi-Co Australia Pty Ltd)</b> 1800 685 899 <a href="http://www.liebherrprofessional.com.au">www.liebherrprofessional.com.au</a>	LKUv1610 LKUv1612 (discontinued) LKUv1613 LKv3910 LKv3912 (discontinued) LKv3913 LKv5710 LKPv1423 LKPv5710 LKPv6253 LKPv6523 FKVSL3612		
<b>MediLine Equipment</b>	125ML	<b>Labcold:</b>	

03 9562 3355 <a href="mailto:info@mediline.com.au">info@mediline.com.au</a>	300ML 400ML ML125GP ML125P ML300GP ML300P ML400GP ML400P ML500GP ML500P	EC75 EC110 EC300 EC400 EC500	
<b>Nuline Refrigeration</b> 03 9562 3355 <a href="mailto:info@nulinerefrigeration.com.au">info@nulinerefrigeration.com.au</a>  <i>*Discontinued – older models</i>	HR200 HR200GP HR200P HR400 HR400G HR400GP HR400P HR600 HR600G HR600GP HR600P HR700GP HR700P	*BEV32P *BEV40P *BEV60P *Mondial Elite BEVPV40 *Elite 310	*GN1420TNP MF1210TNG MF1410TNG *ML400 NLAB2 NLMB2 NLMB/2 P NLM700/1 NLM700-1P NLM1000/P NLM1002 NLMB1006 NM1000/2 *NMM500/1 *NMM1000/2 *TGG-14P *TTGPR14L NLAB215520
<b>Orford Refrigeration Pty Ltd</b> 07 4639 9000 <a href="http://www.orford.com.au">www.orford.com.au</a>	V V D P07A (small) V V D P09A (med) V V D P14A (Large)	FM20S2-R-A FM36S2-R-A FM45S2-R-A	
<b>Quantum Scientific</b>	<b>Sanyo Medi-Cool:</b> MPR-161D(H)* MPR-311D(H)* MPR-514D(H) <i>*Discontinued models</i>	<b>Sanyo Labcool:</b> MPR-721	
<b>Rapini Pty Ltd</b> 03 9555 0848 <a href="http://www.rapini.com.au">www.rapini.com.au</a>	R155-G R300-G WMM100/2P		
<b>Rollex Vaccine Refrigerators</b> Rollex Group (Australia) Pty Ltd 1300 880 441 <a href="http://www.rollexmedical.com.au">www.rollexmedical.com.au</a>  <i>*MATOS-cloud models are a variation of MATOS-eco, which include online monitoring capability and are also compliant.</i>  <i>**Discontinued models</i>	<b>LEC:</b> 350 503 KS290 PESR41 (PE102) PESR47 (PE109) PE202** PESR82 (PE207) PE350 PE357** PE501** PE502** PE507** PE900 PE901** PE902** PE907** PE1202** PE1602** PE1607	<b>Pol Eko–MATOS1*:</b> MATOS-eco 68 R/G MATOS-eco 150 R MATOS-eco 150 R/G MATOS-eco 200 R MATOS-eco 200 R/G MATOS-eco 250 R MATOS-eco 300 MATOS-eco 300 R MATOS-eco 300 R/G MATOS-eco 493 R MATOS-eco 493 R/G MATOS-eco 625 R MATOS-eco 625 R/G MATOS-eco 1460 R/G MATOS-eco PLUS 250 glass/solid door CHL 1+ A CHL 2+ CHL 2+ A	<b>Aria – MATOS:</b> MATOS Aria-Eco Cloud 150 glass/solid door MATOS Aria-Eco Cloud 220 glass/solid door MATOS Aria-Eco Cloud 300 glass/solid door MATOS Aria-Eco Cloud 374 glass/solid door  <b>RK – Koldpia:</b> RK600G RK600S RK1100G RK1100S SGRR-23FD SGRR-23GD SGRR-23SD SGRR-49FD SGRR-49GD

	PEGR41 (PG102) PG202** PEGR82 (PG207) PG357** PG501** PG502** PG507** PG901** PG902** PG907** PG1202** PG1602** PG1607 PEGR158 (PGR151) PEGR273 (PGR273) PGR334** PEGR353 (PGR353) PESR158 (PSR151) PESR273 (PSR273) PESR353 (PSR353) PSR334** PT555	CHL 3+ A CHL 4+ CHL 5+ CHL 5+ A CHL 500 CHL 500 A CHL 700 CHL 700 A CHL 1200 CHL 1450 TROPIC  <b>MATOS Plus Cloud/Eco:</b> 1365R/DT 1365R/GDT 1460R/RT 1460R/GDT 250L Eco Glass 250L Eco Solid 400L Glass 400L Solid	SGRR-49SD  <b>Angelantoni:</b> FRL 360 FRL 360 Glass FRL 500 FRL 500 Glass 700 TN 700 TN-GL 1500 TN 1500 TN-GL  <b>Platinum:</b> LRE 120 LRE 440 PRE 120 PRE 440  <b>Twinbird Portable:</b> SC-BV25
<b>Skope Industries</b>	TME475XL-CK TME100-H TME650-H TME1500-C	B130-PH (discontinued) B170E-PH (discontinued) B310E-PH (discontinued) RCTME1000	LME620/655/190 SK650-2-CK SK1300-2-C SK1500XLP-C
<b>Thermoline Scientific Refrigerators</b> 02 9604 3911 <a href="http://www.thermoline.com.au">www.thermoline.com.au</a>  <i>Note:</i> 1-GD = 1 Glass Door 2-GD = 2 Glass Door 3-GD = 3 Glass Door	TEPR-60-1-GD TEPR-80-1-GD TEPR-140-1-GD* TEPR-145-1-GD TEPR-200-1-GD TEPR-230-1-GD TEPR-250-1-GD TEPR-400-1-GD TEPR-440-1-GD TEPR-500-1-GD TEPR-740-2-GD TEPR-142S-R *Discontinued model	TLMPR140-1-GD TLMPR220-1-GD TLMPR300-1-GD TLMPR350-1-GD TLMPR400-1-GD TLMPR800-2-GD TLMPR900-2-GD TLMPR1215-3-GD	TPR-250-1-GD TPR-520-1-GD TPR-530-1-GD TPR-750-2-GD TPR-900-2-GD TPR-950-2-GD TPR-1086-2-GD
<b>Vaccine Fridges Australia (Frigidum Pty Ltd)</b> 1300 272 464 <a href="http://www.vaccinefridges.com.au">www.vaccinefridges.com.au</a>	Huxford VXH20		
<b>Yili Ramli &amp; Co. Pty Ltd</b> 02 9679 9488 <a href="http://www.vacinerefrigeration.com.au">www.vacinerefrigeration.com.au</a>	VaccineCab YC-280		
<b>Other</b>	<b>Linder and May</b> LMPR-400-1-GD  <b>Tropical Thermaster</b> LG-570GTH (discontinued) TG-570GTH  <b>ARCUS CUI</b>	<b>Baker Refrigeration</b> Custom built  <b>Norcool</b> Super5  <b>RFORD FM45D</b>	<b>Rowson</b> BYY-300  <b>Sanden</b> SPE-0403 SAC-120SAX  <b>MedilineML 300GP</b>